

# **MICROSCOPE ASSESSMENT GUIDELINE FOR ImageStudio SOFTWARE**



# REVISION HISTORY

<b>Manual Version:</b>	A1	<b>Description :</b>
<b>Kit Version:</b>	V1.0 & V2.0	
<b>Date:</b>	May. 2023	<ul style="list-style-type: none"><li>· Initial release</li></ul>

---

<b>Manual Version:</b>	B1	<b>Description :</b>
<b>Kit Version:</b>	V1.0	
<b>Date:</b>	Oct. 2023	<ul style="list-style-type: none"><li>· Added microscope evaluation procedures for Stereo-seq Transcriptomics Multi-immunofluorescence (mIF) solution and Stereo-seq Transcriptomics H&amp;E solution.</li><li>· Added instructions for microscope hardware evaluation.</li><li>· Eliminated the evaluation requirements for blank Stereo-seq Microscope Assessment Chip T.</li><li>· - Modified compatible image formats (only compatible with TIFF/TIF formats).</li></ul>

**Note: Please download the latest version of the manual and use it with the software specific to this manual.**

©2023 Beijing Genomics Institute at Shenzhen.  
All rights reserved.

1. The products shall be for research use only, not for use in diagnostic procedures.
2. The Content on this manual may be protected in whole or in part by applicable intellectual property laws. Beijing Genomics Institute and/or corresponding right subjects own their intellectual property rights according to law, including but not limited to trademark rights, copyrights, etc.
3. Beijing Genomics Institute do not grant or imply the right or license to use any copyrighted content or trademark (registered or unregistered) of us or any third party. Without our written consent, no one shall use, modify, copy, publicly disseminate, change, distribute, or publish the program or Content of this manual without authorization, and shall not use the design or use the design skills to use or take possession of the trademarks, the logo or other proprietary information (including images, text, web design or form) of us or our affiliates.
4. Nothing contained herein is intended to or shall be construed as any warranty, expression or implication of the performance of any products listed or described herein. Any and all warranties applicable to any products listed herein are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. Beijing Genomics Institute, Shenzhen makes no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein.

# TABLE OF CONTENTS

## CHAPTER 1: INTRODUCTION

1.1. Microscope Assessment Channels	3
1.2. Assessment Software and Criteria	3
1.3. Assessment Duration	4

## CHAPTER 2: MICROSCOPE HARDWARE EVALUATION

## CHAPTER 3. MICROSCOPE ASSESSMENT

3.1. Imaging System Recommendations	3
3.2. Imaging Configurations	3
3.3. Input Files and Requirements for Image Registration	4
3.4. Imaging Examples	7

## CHAPTER 4. MICROSCOPE IMAGE ASSESSMENT

4.1. Imaging System Requirements	12
4.2. Pre-assessment Preparation	13
4.3. Microscope Assessment Procedure	14
4.4. Interpretation of Image QC Evaluation Results	31

## CHAPTER 5. INTRODUCTION TO IMAGESTUDIO

5.1. General Information	12
5.2. Installation Instructions	13
5.3. Additional Interpretation of Image Quality Evaluation Results	14
5.4. ImageStudio Software	31



**CRITICAL STEPS:** Pay extra attention for these steps to avoid experimental setbacks or problematic results.



**NOTE:** Additional operation tips and guidance.



**SOLUTION:** Provides a solution to the operation procedure.



**CAUTION:** Proceed with extra care; improper handling or carelessness may cause experimental failure or accidents.

# CHAPTER 1

## INTRODUCTION

STOmics Stereo-seq Solution quantifies mRNA of tissue sections followed by mapping and visualizing the transcriptomic data on an anatomical image. Hence, the tissue image quality determines the success of downstream analysis and spatial clustering of gene expression data. This guidebook provides hardware recommendations, image acquisition and evaluation guidelines, as well as microscope compatibility test procedures. Additionally, a BGI-developed image evaluation and processing software – ImageStudio is also introduced in this guidebook.

## 1.1. Microscope Assessment Channels

- a. Stereo-seq Transcriptomics application(ssDNA-stained nuclei image) -> FITC Channel
- b. Stereo-seq Transcriptomics with mIF co-detection application (DAPI-stained nuclei image & mIF images) -> DAPI Channel
- c. Stereo-seq Transcriptomics H&E application (H&E-stained image) -> Brightfield/Epi-illumination

## 1.2. Assessment Software and Criteria

All microscope images are evaluated for image quality via the "Image Quality Control" module in ImageStudio software, including the scoring evaluation of Trackline clarity and image clarity. All evaluations are based on the standard of Trackline QC score. For detailed scoring rules, please refer to Chapter 4.4 and Chapter 5.3. The image clarity score is for reference only and is not used as a pass/fail standard for image QC.

## 1.3. Assessment Duration

About half a day.

# CHAPTER 2

## MICROSCOPE HARDWARE EVALUATION

## 2.1. Imaging System Recommendations

The microscope hardware evaluation mainly evaluates whether the microscope hardware parameters meet the Stereo-seq experiment requirements. The microscope hardware evaluation contents are as follows:



**It is recommended to make an appointment in advance with your manufacturer's microscope engineer to accompany the evaluation on site.**

Index/ Parameter	Description
XY stage travel distance of microscope	At least 25*75 mm
Focusing approach	Pre-focus map and/or real-time autofocus
Brightfield light	Reflective light source (Epi-illumination)
Objective lens	10X (NA ≥ 0.3)
Fluorescent channel	<ul style="list-style-type: none"> <li>• DAPI (Excitation 358nm, Emission 461nm)</li> <li>• FITC filter cube (Excitation 480/40nm, Emission 525/50nm)</li> <li>• TRITC filter cube (Excitation 545/25 nm, Emission 605/70 nm)</li> <li>• CY5 filter cube (Excitation 620/50 nm, Emission 690/50 nm)</li> </ul>
Camera resolution	High image resolution: ≥ 1800 pixels (height) and ≥2000 pixels (width)
Image bit depth	<ul style="list-style-type: none"> <li>• Grayscale image: 8bit/16bit (fluorescent image)</li> <li>• Color image: 3*8bit/3*16bit (H&amp;E or other brightfield image)</li> </ul>
Exposure time	<ul style="list-style-type: none"> <li>• ssDNA/DAPI: 1 milli sec - 2 sec</li> <li>• H&amp;E: 0.1 milli sec - 100 milli sec</li> </ul>
Background balance	Adjustable background balancing function
Distortion correction	Adjustable distortion correction function
Overlap ratio	Adjustable, no less than 5%
Pixel size	Maximum pixel size ≤ 5 μm/pixel
File format	Capable of viewing or exporting stitched images or Field-of-View (FOV) original images in TIFF format, 8/16 bit
PC requirement	Windows 10 x 64 system, 16G memory or beyond