



**CERTIFICATE OF ANALYSIS**

Nova 8 Calibrator Pack  
Lot Number 729187

Ref Number 15196  
Exp Date 2019-04-18  
Date of Manufacture 2017-10-18  
UDI# (01)00385480151962(11)171018(17)190418(10)729187

Cal A L/N	17291087	Spec	Pass/Fail	Cal B L/N	17298009	Spec	Pass/Fail	Cal C L/N	17298009	Spec	Pass/Fail
Na mmol/L	139.0-141.0		Pass	Na mmol/L	74.5-75.5		Pass	Na mmol/L	139.0-141.0		Pass
K mmol/L	3.97-4.03		Pass	K mmol/L	9.8-10.2		Pass	K mmol/L	3.97-4.03		Pass
iCa mmol/L	0.99-1.01		Pass	iCa mmol/L	1.98-2.02		Pass	iCa mmol/L	0.99-1.01		Pass
iMg mmol/L	0.495-0.505		Pass	iMg mmol/L	0.495-0.505		Pass	iMg mmol/L	1.48-1.52		Pass
pH @ 37C	7.274-7.294		Pass	pH @ 37C	6.730-6.750		Pass	pH @ 37C	7.374-7.394		Pass

  

Reference L/N	17291087	Spec	Pass/Fail
KCl mol/L	1.8-2.2		Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO EN 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Inspector  
Title

11/9/2017  
Date