

ABX Minotrol CRP



LOT

MC 435

CONTROL



(Exp.) **2025-07-05**

(YYYY - MM - DD)

Rev 2

ABX Lysebio

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio																				
			CONTROL					1	TOLERANCES TOLERANCE	CONTROL					2	TOLERANCES TOLERANCE	CONTROL					3	TOLERANCES TOLERANCE
			MICROS			PENTRA		MICROS			PENTRA		MICROS				PENTRA						
			EMI CRP	LC-767G		MS CRP		EMI CRP		LC-767G		MS CRP		EMI CRP	LC-767G			MS CRP					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0		1.9		± 0.4	7.7	7.6		7.9		± 1.0	20.6	20.5		21.3		± 2.2			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.30	2.28		2.24		± 0.15	4.63	4.56		4.66		± 0.18	5.48	5.45		5.53		± 0.25			
HB	HGB	g/dl	6.2	6.2		6.1		± 0.4	13.5	13.5		13.9		± 0.5	17.4	17.8		18.4		± 0.7			
		g/l	62	62		61		± 4	135	135		139		± 5	174	178		184		± 7			
		mmol/l	3.85	3.85		3.79		± 0.25	8.38	8.38		8.63		± 0.31	10.81	11.05		11.43		± 0.44			
HT	HCT	%	16.8	16.9		17.0		± 2.0	37.5	36.9		38.2		± 2.5	48.8	49.1		49.2		± 3.0			
		l/l	0.168	0.169		0.170		± 0.020	0.375	0.369		0.382		± 0.025	0.488	0.491		0.492		± 0.030			
VGM	MCV	µm ³ ; fl	73	74		76.0		± 4	81	81		82.0		± 4	89	90		89.0		± 4			
TGMH	MCH	pg	27.0	27.2		27.2		± 2.4	29.2	29.6		29.8		± 2.4	31.8	32.7		33.3		± 2.8			
		fmol	1.67	1.69		1.69		± 0.15	1.81	1.84		1.85		± 0.15	1.97	2.03		2.07		± 0.18			
CCMH	MCHC	g/dl	36.9	36.7		35.8		± 3.0	36.0	36.5		36.4		± 3.0	35.7	36.3		37.4		± 3.0			
		g/l	369	367		358		± 30	360	365		364		± 30	357	363		374		± 30			
		mmol/l	22.93	22.82		22.25		± 1.86	22.35	22.70		22.59		± 1.86	22.15	22.54		23.22		± 1.86			
IDR	RDW	%	14.0	14.0		15.5		± 3.0	14.0	13.0		13.0		± 3.0	14.0	13.0		11.5		± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	69	70		67		± 20	247	244		243		± 40	452	477		485		± 60			
VPM	MPV	µm ³ ; fl	7.3	7.8		9.0		± 2.0	7.5	7.7		8.3		± 2.0	7.6	7.8		8.3		± 2.0			
LYMPHO	#	1.3	1.3		N/A		± 0.2	2.6	2.5		N/A		± 0.5	3.4	3.3		N/A		± 1.1				
	%	65.0	64.0		N/A		± 8.0	33.5	32.5		N/A		± 6.0	16.5	16.0		N/A		± 5.0				
MONO	#	0.2	0.2		N/A		± 0.2	0.6	0.5		N/A		± 0.5	2.1	1.5		N/A		± 1.4				
	%	9.5	9.5		N/A		± 8.6	7.5	7.0		N/A		± 6.3	10.0	7.5		N/A		± 6.8				
GRANULO	#	0.5	0.5		N/A		± 0.2	4.5	4.6		N/A		± 0.8	15.1	15.7		N/A		± 2.1				
	%	25.5	26.5		N/A		± 9.9	59.0	60.5		N/A		± 9.9	73.5	76.5		N/A		± 9.9				
CRP	mg/dl	0.46	0.51		0.52		± 0.40	2.49	2.64		2.79		± 0.80	4.20	4.23		4.72		± 1.50				
	mg/l	4.6	5.1		5.2		± 4.0	24.9	26.4		27.9		± 8.0	42.0	42.3		47.2		± 15.0				

1300092170
BACK/VERSO
Ref: TEMP-1070 Rev.43