

ABX Minotrol CRP



LOT

MC 429

CONTROL



(Exp.) 2024-01-05

(YYYY - MM - DD)

Rev 1

ABX Lysebio

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio																	
		CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS			PENTRA	MICROS			PENTRA	MICROS				PENTRA					
		EMI CRP	LC-767G		MS CRP			EMI CRP	LC-767G		MS CRP			EMI CRP	LC-767G		MS CRP		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.0		2.1		± 0.4	7.6	7.6		7.9		± 1.0	20.2	20.7		21.1		± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.34		2.31		± 0.15	4.61	4.58		4.61		± 0.18	5.68	5.73		5.76		± 0.25
HB HGB	g/dl	6.1	6.0		5.8		± 0.4	13.5	13.5		13.7		± 0.5	17.8	17.9		18.2		± 0.7
	g/l	61	60		58		± 4	135	135		137		± 5	178	179		182		± 7
	mmol/l	3.79	3.73		3.60		± 0.25	8.38	8.38		8.51		± 0.31	11.05	11.12		11.30		± 0.44
HT HCT	%	16.5	16.4		16.1		± 2.0	37.3	36.6		37.1		± 2.5	49.4	49.3		49.2		± 3.0
	l/l	0.165	0.164		0.161		± 0.020	0.373	0.366		0.371		± 0.025	0.494	0.493		0.492		± 0.030
VGM MCV	µm ³ ; fl	70	70		69.5		± 4	81	80		80.5		± 4	87	86		85.5		± 4
TGMH MCH	pg	26.0	25.6		25.1		± 2.4	29.3	29.5		29.7		± 2.4	31.3	31.2		31.6		± 2.8
	fmol	1.61	1.59		1.56		± 0.15	1.82	1.83		1.85		± 0.15	1.95	1.94		1.96		± 0.18
CCMH MCHC	g/dl	37.1	36.6		36.1		± 3.0	36.2	36.8		36.9		± 3.0	36.0	36.3		37.0		± 3.0
	g/l	371	366		361		± 30	362	368		369		± 30	360	363		370		± 30
	mmol/l	23.03	22.75		22.43		± 1.86	22.45	22.88		22.93		± 1.86	22.37	22.56		22.95		± 1.86
IDR RDW	%	14.5	14.5		16.0		± 3.0	13.0	13.0		13.0		± 3.0	12.0	12.5		12.0		± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	76	80		69		± 20	258	256		254		± 40	495	510		511		± 60
VPM MPV	µm ³ ; fl	8.0	8.7		9.5		± 2.0	8.2	8.4		9.0		± 2.0	8.0	8.2		8.4		± 2.0
LYMPHO	#	1.4	1.3		N/A		± 0.2	2.7	2.6		N/A		± 0.5	3.3	3.1		N/A		± 1.1
	%	65.5	64.5		N/A		± 8.0	35.0	34.0		N/A		± 6.0	16.5	15.0		N/A		± 5.0
MONO	#	0.2	0.2		N/A		± 0.2	0.5	0.6		N/A		± 0.4	1.3	1.4		N/A		± 1.0
	%	7.5	8.5		N/A		± 7.0	7.0	7.5		N/A		± 5.0	6.5	7.0		N/A		± 5.0
GRANULO	#	0.6	0.5		N/A		± 0.2	4.4	4.4		N/A		± 0.5	15.6	16.1		N/A		± 1.5
	%	27.0	27.0		N/A		± 7.0	58.0	58.5		N/A		± 6.0	77.0	78.0		N/A		± 7.0
CRP	mg/dl	0.46	0.49		0.48		± 0.40	2.30	2.41		2.22		± 0.80	4.24	4.35		4.18		± 1.50
	mg/l	4.6	4.9		4.8		± 4.0	23.0	24.1		22.2		± 8.0	42.4	43.5		41.8		± 15.0

Ref: TEMP-1070 Rev.42 1300092170 BACK/VERSO